

MIIG[®]X3

High Strength Injectable
Graft for Compression Fractures

SURGICAL TECHNIQUE



WRIGHT.

MIIG[®] X3

HIGH-STRENGTH INJECTABLE GRAFT FOR
COMPRESSION FRACTURES

surgical techniques as described by
KEITH RASKIN, MD
J. TRACY WATSON, MD



MIIG X3

HIGH STRENGTH
For Compression

The next generation least-invasive grafting option for compression fractures is an injectable calcium sulfate with significantly higher compressive strength. MIIG® 115 (Minimally Invasive Injectable Graft) was the first injectable graft to provide temporary support and resorb at a predictable rate. Now, MIIG®X3 Graft provides even higher compressive strength with the same reliable resorption.

MIIG® X3 Graft

- Greater compressive strength than NORIAN® SRS'
- Injects
- Resorbs
- Remodels into bone




FIGURE 1 |



MIIG® Graft in soft tissues.
MIIG® Graft resorbed.


MIIG® Graft replaced by bone.

MIIG® X3 Graft - Resorbable Salt	VS.	Injectable Calcium Phosphate - Slow Resorbing Ceramic
<ul style="list-style-type: none"> • No Extravasation Issues - Resorbs quickly in soft tissues and synovium 		<ul style="list-style-type: none"> • Ceramic must be removed if leaked into joint space or soft tissues
<ul style="list-style-type: none"> • Faster resorption and replacement by bone (10-12 weeks) 		<ul style="list-style-type: none"> • Resorption rate 2-3 years
<ul style="list-style-type: none"> • Straightforward surgical technique: <ol style="list-style-type: none"> a) Reduce b) Inject graft c) Drill/tap for final hardware 		<ul style="list-style-type: none"> • Backwards technique <ol style="list-style-type: none"> a) Reduce c) Place final hardware b) Inject graft

High-Strength


- MIIG® X3 Graft provides 160% greater compressive strength than MIIG® 115 Graft

FIGURE 2 | Compressive Strength in a Dry Environment



In a dry environment, MIIG® X3 Graft demonstrates compressive strengths comparable to reported bone cement values.

FIGURE 3 | Compression Strength in a Simulated In Vivo Environment



At one hour in a simulated wet environment, MIIG® X3 Graft maintains 160% greater compressive strength than the MIIG® 115 Graft.

Resorbable

- Reliable resorption and replacement by bone

FIGURE 4A |




FIGURE 4B |



A 13x50mm cylindrical defect (outlined) was surgically created in a canine humera model and filled with MIIG® X3 calcium sulfate. Contact radiographs demonstrate an osteoconductive effect through new bone regeneration and subsequent remodeling at 13 weeks.

FIGURE 5 |





Ordering Information

84XS-0405	5cc
84XS-0415	15cc

For Open Reduction, Internal Fixation, Try:



ALLOMATRIX® DR GRAFT

86DR-0300	3cc
-----------	-----



ALLOMATRIX® CUSTOM PUTTY

86XC-0560	5cc
86XC-1000	10cc



LOCON-T® DISTAL RADIUS PLATING SYSTEM

4900-1000

Design Features

- High-Strength
- Injects through small gauge needle
- Hardens in defect
- Replaced by bone
- Temporary support medium
- May be drilled or tapped for hardware

References

- ¹ Reported at 50MPa at 24 hours.
- ² Hand-mixed.
- ³ Vacuum-mixed.
- ⁴ Belkoff SM, et al. The effect of the monomer-to-powder ratio on the material properties of acrylic bone cement. J Biomed Mater Res (Appl Biomater) 2002; 63:396-399. Monomer-to powder ratio was 0.45 ml/g.



Wright Medical Technology, Inc.
 5677 Airline Road
 Arlington, TN 38002
 901.867.9971 phone
 800.238.7188 toll-free
 www.wmt.com

Wright Medical Europe SA
 Commerce Parc
 Gebouw B, 6th Floor Krijgsman II
 1186 DM Amsterdam
 The Netherlands
 Tel +31(0) 205450100

NORIAN® is a registered mark of Synthes Stratec.
 ™Trademarks and ®Registered marks of Wright Medical Technology, Inc.

©2006, Wright Medical Technology, Inc. All Rights Reserved.

SK078-303
 Rev. 3.06