

HIP SYSTEM solutions



CONSERVE[®] Total

HIP SYSTEM with BFH[™] TECHNOLOGY



WRIGHT.


CONSERVE® Total

HIP SYSTEM with BFH™ TECHNOLOGY

Designed to mimic the anatomic size of the femoral head, the CONSERVE® Total Hip System with BFH™ technology provides a low wear, highly stable, versatile modular system to better reproduce human kinematics in both primary and revision surgical settings.


With or Without HA

Shells range in size from 42-64 in 2mm increments

147 to 169 degrees ROM with a jumping distance of 16-25mm

Shells are designed to provide a 6mm difference between the reamed acetabular and the femoral head component

Femoral components range in size from 36-56 in 2mm increments



ACETABULAR CUPS

| CATALOG NO. HA | CATALOG NO. Non HA | CUP OD | CUP ID |
|-------------------|-----------------------|--------|--------|
| 38HA-3642 | 3802-3642 | 42 | 36 |
| 38HA-3844 | 3802-3844 | 44 | 38 |
| 38HA-4046 | 3802-4046 | 46 | 40 |
| 38HA-4248 | 3802-4248 | 48 | 42 |
| 38HA-4450 | 3802-4450 | 50 | 44 |
| 38HA-4652 | 3802-4652 | 52 | 46 |
| 38HA-4854 | 3802-4854 | 54 | 48 |
| 38HA-5056 | 3802-5056 | 56 | 50 |
| 38HA-5258 | 3802-5258 | 58 | 52 |
| 38HA-5460 | 3802-5460 | 60 | 54 |
| 38HA-5662 | 3802-5662 | 62 | 56 |
| 38HA-5664 | 3802-5664 | 64 | 56 |

FEMORAL HEADS

| HEAD OD | CATALOG NO. -3.5 mm | CATALOG NO. 0 mm | CATALOG NO. + 3.5 mm |
|---------|------------------------|---------------------|-------------------------|
| 36 | 2600-0025 | 2600-0026 | 2600-0027 |
| 38 | 3801-3804 | 3801-3800 | 3801-3835 |
| 40 | 3801-4004 | 3801-4000 | 3801-4035 |
| 42 | 3801-4204 | 3801-4200 | 3801-4235 |
| 44 | 3801-4404 | 3801-4400 | 3801-4435 |
| 46 | 3801-4604 | 3801-4600 | 3801-4635 |
| 48 | 3801-4804 | 3801-4800 | 3801-4835 |
| 50 | 3801-5004 | 3801-5000 | 3801-5035 |
| 52 | 3801-5204 | 3801-5200 | 3801-5235 |
| 54 | 3801-5404 | 3801-5400 | 3801-5435 |
| 56 | 3801-5604 | 3801-5600 | 3801-5635 |
| 56 | 3801-5604 | 3801-5600 | 3801-5635 |

MATCH with the following Wright Stems:

Any PROFEMUR® Series Stem

ANCA-FIT®

PERFECTA®

HELIANTHE®



Wright Medical Technology, Inc.

5677 Airline Road
Arlington, TN 38002
901.867.9971 phone
800.238.7188 toll-free
www.wmt.com

Wright Cremascoli Ortho SA

Zone Industrielle la Farlecle
Rue Pasteur BP 222
83089 Toulon Cedex 09
France
+ 33.49.408.7788 phone