

extremities solutions



WRIGHT.

E-CENTRIX[®]
Ulnar Head

Anatomically Designed. For Replacement of the Distal Ulnar Head.

E-CENTRIX™ ULNAR HEAD

INDICATIONS

Replacement of the distal ulnar head for disorders of the distal radioulnar joint in rheumatoid, degenerative and post-traumatic arthritis presenting with:

- Pain and Weakness of the Wrist Joint (not improved by nonoperative treatment);
- Instability of the Ulnar Head with Radiographic Evidence of Dislocation or Erosive Changes of the Distal Radioulnar Joint;
- Failed Ulnar Head Resection (e.g. Darrach Resection);
- Primary Replacement After Fracture of the Ulnar Head or Neck;
- Revision Following Failed Ulnar Head Arthroplasty

Suture Holes allow for anatomic soft tissue attachment including triangular fibrocartilage complex

Cement or cementless options

Modular design for optional sizing including revision length resection

Eccentric anatomical head articulation

Plasma porous ingrowth surface for soft tissue integration stability

BENEFITS

- Pain Relief of Radioulnar Impingement and Radioulnar Instability
- Restores Normal Wrist Appearance
- Eccentric Anatomical Design of the Head Avoids Overstuffing the Joint and Restores Kinematics
- Superior Stable Suture Attachment Options for TFCC and Joint Capsule
- More Stem Options For Better Press-Fit and Revision Applications
- Head Made of Cobalt Chrome and Stem of Titanium

ORDERING INFORMATION

Ulnar Heads

596H1615	16 x 1.5mm Offset
596H1630	16 x 3.0mm Offset
596H1815	18 x 1.5mm Offset
596H1830	18 x 3.0mm Offset
596H2015	20 x 1.5mm Offset
596H2030	20 x 3.0mm Offset

Ulnar Stems

596S0055	5.5mm Standard
596S0065	6.5mm Standard
596S0075	7.5mm Standard
596S0085	8.5mm Standard
596S0100	Ulnar Stem Extender



Wright Medical Technology, Inc.

5677 Airline Road
Arlington, TN 38002
901.867.9971 phone
800.238.7188 toll-free
www.wmt.com

Wright Cremascoli Ortho SA

Zone Industrielle la Farlede
Rue Pasteur BP 222
83089 Toulon Cedex 09
France
011.33.49.408.7788 phone