



WRIGHT.

HUNTER™

Active and Passive Tendon Implants

A consistently reliable method of salvaging scarred tendon systems

HUNTER™

Active and Passive Tendon Implants

A Consistently Reliable Method of Salvaging Scarred Tendon Systems

Designed for temporary use in reconstruction of extension and flexor tendon based on a two stage procedure.

Active tendons allow use of hand during rehabilitation thereby:

- Increasing patient motivation to comply with therapy.
- Reducing loss of range of motion before autogenous graft implantation.
- Reducing proximal juncture adhesion at the site of motor tendon sheath.
- Allowing active evaluation of the pulley system before autogenous graft placement.

Bone Screw Plate available providing firm fixation of distal attachment further reducing the risk of detachment and migration

Manufactured with polyester reinforcement providing added strength and reducing risk of detachment and migration

Rods create a pseudosynovial sheath which later lubricates and nourishes an autogenous graft.

HUNTER ACTIVE TENDON IMPLANT (AT00-KIT B)

4 MM WIDTH



16 cm Long	AT164000
18 cm Long	AT184000
20 cm Long	AT204000
22 cm Long	AT224000

4 MM WIDTH



16 cm Long	ATBC1640
18 cm Long	ATBC1840
20 cm Long	ATBC2040
22 cm Long	ATBC2240
26 cm Long	ATBC2640

4 MM WIDTH



16 cm Long	ATPC1640
18 cm Long	ATPC1840
20 cm Long	ATPC2040
22 cm Long	ATPC2240
26 cm Long	ATPC2640

4 MM WIDTH



16 cm Long	ATDC1640
18 cm Long	ATDC1840
20 cm Long	ATDC2040
22 cm Long	ATDC2240
26 cm Long	ATDC2640

HUNTER PASSIVE TENDON IMPLANT (PTTR-KIT B)

25 CM LENGTH



3 mm Wide	PT300000
4 mm Wide	PT400000
5 mm Wide	PT500000
6 mm Wide	PT600000

24.5 CM TENDON ROD LENGTH



2 mm Wide	TR200000
3 mm Wide	TR300000
4 mm Wide	TR400000
5 mm Wide	TR500000
6 mm Wide	TR600000

WRIGHT.

Wright Medical Technology, Inc.

5677 Airline Road
Arlington, TN 38002
901.867.9971 phone
800.238.7188 toll-free
www.wmt.com

Wright Cremascoli Ortho SA

Zone Industrielle la Farlecle
Rue Pasteur BP 222
83089 Toulon Cedex 09
France
011.33.49.408.7788 phone

™Trademarks and ®Registered marks of Wright Medical Technology, Inc.

©2005 Wright Medical Technology, Inc. All Rights Reserved.

SO 773-798
Rev. 5-05