



RAYHACK® Kienbock Radial Shortening System Saw Blade Summary



Two individually packaged sterile saw blades are provided with each RAYHACK® Ulnar and Kienbock system. Both saw blades 40100110 and 110 are compatible with the LINVATEC/HALL®/MICROAIRE® saws. Part numbers 40100410 and 410 saw blades are compatible with STRYKER® saws.

The provided saw blades are specifically designed for use with the RAYHACK® Ulnar and Kienbock systems. Using a generic blade is considered off-label and the responsibility of the user.

The below chart is for reference only. Due to manufacturer's specification changes Wright cannot guarantee the saw blades will fit all systems.

		WRIGHT PART NO	PNEUMATIC	ELECTRIC	BATTERY
			1000-100/1950	1000ET/1950	6642B/6672
			1000-100/1955	000ET/1955	6642B/6670 (Battery SmartDriver)
				1641/6672 (Electric SmartDriver)	
 MICROAIRE® and HALL® (LINVATEC) compatible	ZS110*	40100110 110	5053-011 MICRO 100™	6020-022 MICROPOWER® (Pencil Style)	Pro6200/Pro2043 MPOWER™
			Pro6150/Pro2043 PNEUMATIC POWERPRO®	6021-022 MICROPOWER® (Pencil Style)	Pro6200/Pro2043 MPOWER™
				5020-022 MICROCHOICE® (Pencil Style)	Pro5100/Pro2043 POWERPRO®
				Pro6100/Pro2043 (Electric POWERPRO®)	Pro5100/Pro2043 POWERPRO®
 STRYKER® compatible	SP440 R*	40100410 410		Command II 2296-34	Cordless Driver 4200 with 4100-400 Saw
				TPS 5100-34	
				CORE 5400-34	
				5400-99/4100-400 (Core Driver with Saw)	

NOTE: Right angle grip of battery drivers is more difficult to use compared to pencil-type grips.

*MICROAIRE® Part No.